


The best choice in beauty equipment

GENERAL CATALOGUE






Experts in beauty technology since 1978


RÖS'S company is expert in manufacturing beauty technology since 1978. Our core values are high quality, constant technology innovation and excellence in aftersales and customer services. We have a consolidated experience in OBM & OEM manufacture.

Since its foundation, by Carlos Moliner, Vicenç Gil and Francesc Esteve, RÖS'S has kept alive its commitment: to give an all-round service and to provide beautician professionals with effective technological solutions. The company is today consolidated as an absolute reference for technological innovation in the sector.

The concern of our team of professionals, made up of beauticians, engineers, physiotherapists and beautician doctors, is to seek new safe and efficient techniques to offer the right solutions for each business.

Thanks to our aspiration for international growth, we now have a presence in more than 30 countries worldwide. For all our present and future clients, we place our facilities at your disposal, with over 5,500 square metres devoted to research, design, production, advice and training in the latest technological advances.





Our values. Our company

Excellence in manufacturing

As a company with ISO certification as "Designers, manufacturers and suppliers of electromedical and esthetic electrotherapeutic equipments", we pay great attention to quality protocols during the development and manufacturing process, complying with the ISO 9001:2008 international standard.

Passion for innovation

RÖS'S aims for new combinations of energy and intuitive applications to optimise results. Our constant concern is to research to find new solutions that are as easy as they are surprising.

Always close beside the professional

Our team works with a vocation to provide service, close at hand to provide advice at all times when it comes to taking decisions. Sharing our knowledge and training professionals are the main options we see for the future of the company. We accompany our clients day by day in their growth. Our customer service telephone number, (34) 934 266 017 is at your disposal.

Beauty and wellbeing

It is our responsibility to work with rigour, always according priority to the safety of all our treatments. That is why we place our trust in a technological line controlled by qualified engineers, with non-invasive application protocols developed under medical control.

A world of possibilities in aesthetic technology



The RÖS'S model is based on the development and creation of a balanced range of designer products for beauty salons, offering intelligent use of energy and effective treatment.


The success of RÖS'S technological solutions lies both in the company's ability to recognise and adapt to changing trends, as well as in the development of innovative products. Our commitment is to develop reliable, innovative, effective and easy-to-use products that offer a unique method of treatment.


The whole range aims to meet the basic specifications which make our products reliable and exclusive to our clients.

- 1. EU MANUFACTURING
- 2. EXCLUSIVE TECHNIQUES AND APPLICATIONS
- 3. RELIABLE TECHNOLOGY
- 4. PROVEN AND DOCUMENTED RESULTS

Choosing professionalism and service is opting for differentiation and success.









THE PHYSIO-ACTIVE MASSAGE SYSTEM 100% NATURAL

Rollaction

Rollaction System is the revolutionary and exclusive physio-active body massage system inspired by the masseur's professional movements. Rollaction provides strength and energy for the most in-demand body treatments, increasing efficacy with minimum effort.

Manual care and technology have been combined together to provide a versatile tool with proven results.

Our work system kick-starts the metabolism without being invasive, respecting the nature of the vascular, muscle and skin systems. It restores the balance and harmony of the body.

How does Rollaction work?

Rollaction works through its roll-heads which perform a series of high-penetration rotations. This rotation causes friction in the form of kneading on the skin and tissues which, when combined with pressure, simulates the effect of a professional's hands and maximises their effect. The controller also projects **stimulating infrared light** to improve results.



TREATMENTS

- ANTI-CELLULITE
- SCULPTING
- TONING
- CIRCULATION-STIMULATING
- FAT-REDUCING.
- DECONTRACTING
- SPORTS MASSAGE





- Firm skin: Unlike other systems,
 Rollaction acts on deep-rooted cellulite
 and removes unsightly orange-peel skin.
- Sculpting: The rotating action dissolves accumulated fat and facilitates its drainage via the lymphatic channels. With Rollaction, getting back in shape becomes a reality.
- **Toning:** Muscle spasms compress the capillaries nourishing the tissue and stop the blood from reaching it correctly. Through deep physio-active massage, the additional supply of oxygen and nutrients nourishes and restores the muscle fibres.
- Vascular system: Thanks to the mechanical stimulus, the muscle fibres in the blood vessel walls are activated. The improvement in the circulatory system due to the rise in blood supply increases metabolic activity.
- Wellness effects: All of our tissues (muscle, bone, connective, etc.) require a good blood supply. Incorrect postures, a sedentary lifestyle, muscle lesions, etc. cause our muscles to spasm, obstructing the paths nourishing oxygen and causing significant imbalances in the locomotor system and even organ dysfunction.

Rollaction tacts on entire muscle areas and, thanks to its rotating movement and gradual pressure, it restores circulation and nerve impulses in the affected areas.








The perfect ally for the masseur and physiotherapist

Rollaction

SUCTION


Traditional systems work by clamping or sucking the skin, causing direct mechanical traction on the surface tissues, without reaching the deep layers.





COMPRESSION

Rollaction Rollaction acts on all deep tissues without causing any traction, but rather through regular and controlled compression on all tissues, from the skin to the deep muscles.



Two levels of penetration Maximum efficacy

AT THE DERMAL LEVEL

- | Anti-cellulite programme
- | Fat-reducing programme
- | Toning/sculpting programme
- | Circulation-stimulating programme

AT THE MUSCLE LEVEL

- | Decontracting programme
- | Sports massage programme

STUDY 1 ROLL ACTION

RÖS'S has conducted a thorough in-booth study of the Rollaction system. The study has been conducted in a sample of 50 women with beauty problems caused by circulatory deficiencies.

Mean reduction in volume in cm					
Waist	Hips	Buttocks	Thighs	Knees	
8,16	8,58	6,82	5,44	1,3	

A reduction in volume of up to three dress sizes is seen with reduction of up to 12 cm on the waist and hips.

The texture of the skin is substantially improved thanks to the elimination of orange-peel skin and to tissue firming.

TRIPLE ROLLACTION RESULTS

- | Vascular and lymphatic action: increased arterial blood supply and venous and lymphatic drainage
- I Skin action: unblocking of the vessels and increased blood supply to all tissues
- | Muscle action: it increases the supply of nutrients and oxygen, toning the muscle fibres and eliminating spasms





Press G3

PRESSURE-THERAPY, THE GREAT ALLY OF BEAUTY, WELLBEING AND HEALTH

Pressure therapy is an essential technique to help professionals from beauty and health care areas. Thanks to its proven efficacy, pneumatic massage has become one of the most popular treatments.

Lymph drainage, better circulation and relaxation are just some of the benefits.

RÖS'S has designed 9 pre-established programs, 6 for the legs and 3 for the arms. 9 free programs for customisation.

How does Press G3 work?

Bad habits and today's lifestyle carry the risk of physical and nervous disorders and anxiety. Presso-therapy is a pneumatic massage system based on controlled compression and decompression which prevents the development of circulatory problems and increases relaxation. Medically proven effectiveness.

Application and benefits

- | Stimulation of venous return
- | Prevents the development of varicose veins
- I Ideal for post-pregnancy care
- | Detoxification after surgery (with a doctor's prescription)
- | Prevents the development of cellulite.

General characteristics

| Superimposed air chambers for better drainage

Intensity regulator and pressure time individualised by section

Intelligent pressure test

| Termination of treatments by cycle or time.



PROGRAMS FOR LOWER **EXTREMITIES:**

- LEG DRAINAGE
- LEG ACTIVATOR
- ORANGE PEEL SKIN
- POST-SLIMMING
- TIRED LEGS
- MICROACTIVATOR

PROGRAMS FOR UPPER **EXTREMITIES:**

 POST-MASTECTOMY (only with medical prescription)

15

- ARMS ACTIVATOR
- ARMS DRAINAGE





LATEST-GENERATION THERMO-STIMULATION

TEI System

Infrared band and variable current equipment for joint heat and muscle stimulation treatment.

RÖS'S pioneered the development of thermo-stimulation and its team of engineers, in collaboration with beauty professionals, has designed the most sophisticated application bands to provide customers with the best possible results and maximum safety.

The pleasure of technology Reduction of volume in just four sessions

Two basic professional beauty-treatment agents, heat and stimulation, are the key factors in TEI System. Thanks to deep infrared heat, the adipose cells in fat tissue are heated, aiding their liquefaction by electric current. The controlled emission of heat combined with different types of electric current guarantee the results.

18 electro-currents combined in 21 body programmes



Highly adaptable strips

TREATMENTS

- SOFT CELLULITE
- HARD CELLULITE
- OEDEMATOUS CELLULITE
- FIRMING
- INTENSE FIRMING
- POST-PREGNANCY
- EXCESS WEIGHT
- LOCALISED EXCESS WEIGHT
- CIRCULATORY
- DRAINING
- REDUCTION
- SCULPTING
- TONING
- RELAXATION
- TENS
- LIPOLYTIC
- PASSIVE GYMNASTICS
- CALMING PASSIVE GYMNASTICS
- SEQUENTIAL BOOSTING
- LIPO-DRAINING
- UPPER TRUNK





Ergonomic

cables-carrier

TEI System

Advantages of using two techniques simultaneously

- | The feeling of electro-stimulation is more pleasant thanks to the sensation brought by thermotherapy
- | Shorter application times
- Increased electrical conductivity from sweating
- | Benefits of both static and dynamic heat

Benefits derived from TEI System

- | Lipolytic effect from penetration of infrared heat from 1 cm to 3.5 cm
- | Stimulation of blood and lymphatic circulation
- | Increased metabolic activity
- Improved vascular system and venous return
- Increased cell tropism, healing capacity and connective tissue elasticity
- | Analgesia and muscle relaxation
- | Muscle firming with all treatments

Latest-generation strips

The easily applied thermo-stimulation strips are the key to **TEI System** treatment. Flexible, adaptable, tough and safe, their laminar circuit provides completely uniform infrared emission, while their flexible cover adapts to the body to ensure full contact. .

Supplied with 12 long-lasting strips which adapt to the body. Rapid application. Easy to clean. Four different sizes. There is also a Courtesy Kit to cope with any incident.



Easy application of conductive gel

Thermo-stimulation strips designed to wrap the entire muscle group

Silicone coating. Very flexible and biocompatible to avoid allergic reactions



THERMO-STIMULATION, TWO TREATMENTS IN ONE

Kestos T2

The joint effect of infrared heat therapy and electro-stimulation enables the simultaneous application of two treatments in one. This technique, created by RÖS'S, provides proven results in record time.

Thermo-stimulation is one of the star treatments in professional beauty therapy. **Kestos T2** provides you with unique freedom, as the bands work by themselves while you make the most out of your time.

Kestos T2 is an excellent ally in body treatments; solutions for your everyday needs and proven cost-effectiveness.

How does Kestos T2 work?

Thermo-stimulation, a classic beauty treatment, makes the most of heat and electro-stimulation. Infrared heat penetrates as far as the adipose tissue. Sweating facilitates the conduction of electric currents which have a percussion effect, promoting its breakdown and drainage, while firming the tissues.

Thermo-stimulation is a method which makes the most of heat and electro-stimulation. The effects of both techniques are enhanced:

- I the heat generates sedation, making electro-stimulation even more pleasant while aiding the breakdown of cellulite
- I the stimulation firms the tissues, preventing loose skin in volume reduction or anti-cellulite treatments

Flexible, strong and


easy-to clean latest-

generation strips

TREATMENTS

- SOFT CELLULITE
- HARD CELLULITE
- OEDEMATOUS CELLULITE
- FIRMING
- INTENSE FIRMING
- POST-PREGNANCY
- OVER-WEIGHT
- LOCALISED EXCESS WEIGHT
- CIRCULATORY
- DRAINING
- REDUCTION
- SCULPTING
- TONING
- ANALGESIC EFFECT





US₁

INTENSE ANTI CELLULITE ACTION **TECHNOLOGY**

RÖS'S has designed US1, the ultrasounds equipment, specific for anti cellulite, reduction and drainage body treatments. Its emission frequency, of 1MHz, penetrates between 4-5cm. It's a powerful allay for body treatments, like excess volume or anti aesthetic cellulite.

US1 action has three synergic effects

- | Thermal action: The local increase in tissue temperature activates the metabolic processes, which result in an additional vitality in the tissues.
- | Mechanical action: The mechanical effect of "percussion" decompresses the adipose tissue.
- | Sonophoresis: Very efficient system for the transdermal transport of active substances without them having to be ionisable.

Adjustable parameters to customise the treatment according to the treatment:

| Emission time.

I Three power levels.

| Continuous or pulsed. Three pulsed



The 60mm diameter ultra sound waves allow work to be done in large areas in a quicker, easier and more homogeneous way.

TREATMENTS

- BODY SCULPTING
- REDUCING AND PREVENTING CELLULITE
- CIRCULATION IMPROVEMENT
- NATURAL ELIMINATION OF RETAINED LIQUIDS
- SKIN ASPECT IS GENERALLY IMPROVED





BEAUTY SOLUTIONS BEYOND DIATHERMY

Binary Premium

Binary Premium is the multifunction device for aesthetic facial and body treatments. Its mthod is based on the principles of complementarity and the synergy of cell regenerative diathermy with electro-stimulation.

DIATERMIA

Binary Premium emits at a frequency of 448 kHz, combining capacitive and resistive applications. Its action re-balances the electro-cellular potential in all tissues, triggering biological changes in the body at the biochemical level. The effects of arterial vasodilation enhance the permeability of the cell membrane, increasing oxygenation and nutrition and increasing venous drainage, which means decreased oedema in inflammatory processes. As a result of this process, deep tissue regeneration is initiated, increasing the activation of neocollagenesis and the renewal of dermis, muscle, cartilage and bone.



PRE-ESTABLISHED PROGRAMS

- FIRMING
- REVITALISING
- UPPER INTENSIVE
- MIDDLE INTENSIVE
- LOWER INTENSIVE

BODY TREATMENTS

- ANTI-CELLULITE PROGRAMS
- REDUCTION
- FIRMING

FREE PROGRAMMING FACIAL TREATMENTS

- EYE PUFFINESS
- DOUBLE CHIN
- EXPRESSION LINES
- FACIAL CONTOUR RECOVERY
- POST FACIAL SURGERY

FACIAL TREATMENTS BODY TREATMENTS

- LYMPHATIC DRAINAGE
- BREASTS
- STRETCH MARKS
- MUSCLE CONTRACTURES
- SCAR RECOVERY
- POST BODY SURGERY





DIATHERMY

I CAPACITIVE application

During the application of diathermy via capacitive electrodes, the energy is transmitted as a displacement current by an electrode dressed in an insulating sleeve. The bipolar molecules, in particular water, vibrate and change direction at the same frequency as the current to which they are subjected. As such, the structures most strongly affected are those that are the most vascularised, or those with a higher water content: dermis and muscle.



I RESISTIVE application

During the application of diathermy via resistive electrodes, the energy is transmitted as a conduction current by a metal electrode in direct contact with the skin.

This causes an intense internal temperature increase, which will be better for those structures that are more resistant to the passage of the electric current: fat, dense tissues and fibrotic tissues.





Binary Premium

ELECTRO-STIMULATION

I ELECTROLIFTING

A combination of specific micro-currents that stimulate the skin and muscle microstructures.

- Activation of microcirculation and general metabolism with a revitalising effect for tissues.
- Recovery of muscle tone: relaxation of contracted muscles.



Specific low-frequency currents that act on the adipocytes and the fats they contain, causing their destruction and facilitating their removal. At the same time, it activates the circulatory system to facilitate the removal of fat content.

The hands-free electro-stimulation bands synchronize two processes: the degradation of fat and its drainage.

I ELECTRO-FIRMING

Excitomotor currents that trigger muscle contractions, indicated for muscle rehabilitation and recovery treatments and to stimulate muscle firmness in muscles that have lost their tone.





Pointer-electrodes facilitate precision in each application manoeuvre.



Capacitive angled handle + 3 electrodes



Resistive angled handle + 3 electrodes

SOFTWARE EXPANSION AND ACCESSORIES RANGE

O i i

Pelvic floor kit

Hand-free resistive bracelets



Double Bracelets for Systemic treatments







TERMOBEAUTY IS DIATHERMY EMPLOYING 448 KHZ

TermoBeauty

Any unsightly blemish or disease results in an alteration of cellular electrical potential. The consequences of this affect our morphology very differently, so we need a precise energy that activates or slows down the cell pulse. This global rebalancing and regenerating action progressively recovers damaged tissues in the body. This is science applied to regeneration and recovery.

The 448 kHz high-frequency emission signal is the most appropriate and scientifically endorsed signal to achieve intercellular ion mobilisation and action from the heart of the cell. Its two modes of application, capacitive monopole and resistive monopole, allow us to reach all tissues, from the dermis to the bone tissue.

General treatments

Round electrode



Resistive treatments

Capacitive treatments

Precise treatments Concave electrode



- EYE BAGS
- DOUBLE CHIN
- COUPEROSIS
- ROSACEA
- SEBORRHEIC REGULATION

BODY SOLUTIONS

- BREAST FIRMING
- STRETCH MARKS IN FORMATION
- SCAR RECOVERY
- MARKED LAXITY
 (ABDOMEN, ARMS)
- KELOIDS
- WELLNESS RECOVERY
- CAPILLARY THERAPY

PHASE 1	Capacitive biostimulation	More specific action in dermis and muscle. Induction of local hyperthermia activates collagen and elastin generation.
PHASE 2	Resistive biostimulation	More specific action in dense, fatty and fibrotic tissues. Optimal for adipose tissue reduction treatments.





FOR FIRMING THE SKIN AND DECONGESTION OF TISSUES

Mesobiolift RF

Dermal treatments were until now programmed in phases: firstly, stimulating vasodilation and then causing nutrients to penetrate down inside the skin.

With Mesobiolift RF, we achieve both objectives at the same time. This new system acts in a synchronised way from the surface and from the interior of the skin, with the objective of actuating the mechanisms that boost cellular regeneration.

How does Mesobiolift RF work?

The endothermal energy generated by Mesobiolift RF results from a fusion of electromagnetic energy and direct current.

This new emission of energy is a powerful generator of controlled internal heat with the ability to transdermally introduce active ingredients into the skin.

Professional protocol

PHASE | Scanning from 1MHz

Cutaneous

CControls the progressive increase in skin temperature. Keeps all cutaneous layers stimulated and

PHASE II Endothermal wave

Cellular activator

Treatment in synergy of RF and active ingredients penetration

PHASE III Radiofrequency at 3MHz

Efficacy booster

Circulatory activator to boost takeup of active ingredients




FACIAL SOLUTIONS

- SKIN FIRMING
- FACE OVAL CORRECTION
- PENETRATION OF ACTIVE INGREDIENTS
- CUTANEOUS RESTORING

BODY SOLUTIONS

- SKIN TIGHTENING
- TISSUES DECONGESTION
- PENETRATION OF ACTIVE INGREDIENTS





Nasolabial folds

Facial treatments

In facial treatments, the resistance offered by the tissues to the passing of current leads to a controlled increase in the tissue temperature in the dermis, while leaving the epidermis intact. This endothermal increase has a dual consequence:

Immediate activation of the blood circulation

Reception of Supreme Cell active ingredients and activation of cellular neogenesis

Contraction of the collagen fibres

Immediate effect: tensing of collagen fibres in the subcutaneous tissue

Progressive effect: formation of new elastic fibres

IMPROVEMENT OF SKIN BIO-REDENSIFIYING NET

Firming, face oval correction, improvement of general elastic tone, improvement of expression lines







FACIAL AND BODY TREATMENTS

Mesobiolift RF



RÖS'S conductive emulsions for endothermal wave

The endothermal wave of itself stimulates the growth of the stem cells, but in addition the ionic interchange that it unleashes enhances the permeability of the cellular membrane, which means that penetration of active ingredients is more much easier.

RÖS'S introduces different options of conductive emulsions for endothermal wave enriched with plant active ingredients encapsulated in liposomes.



RF Firming Emulsion

Emulsion enriched with equisetum, hop extract, ivy, hammamelis and collagen, a plant compound for effectively treating loss of skin elasticity. It stimulates tissue drainage and has firming properties on the connective tissue.

IM Decongestant Emulsion

Emulsion enriched with menthol extract, taken from species of mint with anti-inflammatory, antiseptic and analgesic properties, and calmosensine, a peptide which reduces the perception of unpleasant sensations.

Supreme CELL

Emulsion enriched with plant stem cells extracted from the nucleus of the argan fruit, a recognised cellular regenerator. It activates neogenesis of the skin cells, reconstructing the net that holds our cutaneous tissue in place. Optimum for skins devitalised through lack of nutrients and now beginning to show signs of ageing.

Body treatments

The activator power of radiofrequency induces a proliferation of dermal fibroblasts and improves the arteriovenous system of the tissues. It's one of the most suitable energy forms for treating flaccidity in difficult zones such as the inner side of the thighs and on the abdomen or arms.

With the Mesobiolift system we apply a 1Mhz radio frequency and virtual mesotherapy in a single movement, causing plant active ingredients to penetrate down inside the skin and enhancing the treatment results.

FIRMING TREATMENT.

Skin revitalization and regeneration.

DECONGESTANT TREATMENT.

Stimulation of blood circulation on the area.




Buttocks



Abdomen



Arms



ANTI-AGEING CARE WITH RADIO FREQUENCY

My Sekret

My Sekret emits a radio frequency at 3MHz, which is a powerful heat generator with applications for cutaneous firming of the face and body.

How does My Sekret work?

The aesthetic radio frequency consists of applying electromagnetic energy to provoke controlled heating of the different skin (dermis) layers. This controlled temperature increase is produced by tissue resistance to the current. Its application helps the formation of new collagen and elastin, as well as an improvement in circulation in the treated area.

Facial Application

As time goes by, tissue cell renewal becomes slower and the skin loses elasticity because the fibres responsible for recovery, such as collagen and elastin, among others, weaken. The skin gets thinner and lines tend to appear. The latest radio frequency technology aims to prevent and recover skin smoothness and firmness, increasing the cell metabolism and regenerating new elastic fibres.

- Improvement of the dermal elastic tonicity.
- | Stimulation of collagen production.

Body application

Various endocrine, genetic, hormonal or neurovegetative factors favour the appearance of conditions like flaccidity. With radio frequency we detect a gradual improvement in skin quality, as it becomes firmer and smoother.

- Supply of nutrients and oxygen.
- It improves lymphatic draining for the elimination of waste substances.

body and facial appliances.

FACIAL TREATMENTS

- RENEWAL OF SKIN TISSUE
- DERMAL FIRMING

BODY TREATMENTS

• CUTANEOUS FIRMING

Coaxial bipolar resistive

interchangeable electrode for

applicator with an

TISSUE DECONGESTION

Customisable parameters on a touch screen

- | Time of application.
- Application head: body or facial.



PROFESSIONAL CARE FOR THE SKIN

Complet

Complet is the perfect equipment for professional basic facial treatments. Its techniques can be combined to obtain multiple treatment options.

Whatever the skin type, **Complet** has the right solution for facial cleansing and treatment.

Never have there been so many options in such a small space. Useful, practical and versatile.

How does Complet work?

- 1 Suction: Suction pad. Clean skin. Two types of pad, round and flat, tailored to your needs.
- I Deep cleansing: Deep cleansing by galvanic current for treatments on greasy or oily skins. It regulates sebaceous secretion. Vibrant, oxygenated skin.
- I Facial revitalisation by means of micro-stimulation: Emission of micro-currents which act internally to boost tissue oxygenation.
- I Spray: Micro-spray of treatment tonics. Moisturised and toned skin
- I Direct or indirect high frequency: Relaxing or mechanical stimulating. Product penetration thanks to thermal effect. Clean, decongested skin. Four different electrodes.
- Brossage: Brush to remove accumulated impurities and oil with a soapy lotion. Clean, well cared for skin. Four types of brush.
- I lonisation: Introduction of ionisable cosmetic substances through the skin. Moisturised and well cared for skin. Two roller-electrodes.
- I **Electrolifting:** Emission of micro-currents which act on muscle myofibrils, preventing the development of loose skin and expression lines.

TREATMENTS

- SUCTION
- SPRAYING
- BROSSAGE
- IONISATION
- DIRECT HIGH FREQUENCY
- INDIRECT HIGH FREQUENCY
- ELECTRO-LIFTING
- MICRO-STIMULATION
- IONISATION+LIFTING

Spraying, high frequency, ionisation and electro-lifting





EPILATION

Depilplus 3

TREATMENTS

 PERMANENT REMOVAL OF UNWANTED HAIR



DEFINITE REMOVAL OF ALL UNWANTED BODY HAIR

Electric hair removal is the only technique that guarantees the elimination of all non-desired hair regardless of its colour, thickness or location. This technology works with controlled emissions of galvanic and high frequency currents. It is recommended for small areas, residual hair after photoepilation or in cases where traditional hair removal is not possible.

DEPILPLUS 3 is the electro-aesthetic equipment, which includes the three most effective techniques in electric hair removal: thermolysis, electrolysis and the combination of both (system blend).

Parameters that may be personalised and which allow us to adapt the application according to the type of body hair.

- | Technique.
- I Time.
- | Power degree.

Technique	Emission	Action	Choice of technique
Electrolysis	Galvanic current	Destruction of the dermal papillae by (alkaline) chemically burning. A return electrode is required.	A high level of efficiency and a medium level of inflammation after treatment.
Thermolysis	High frequency	Destruction of the dermal papillae by thermal action.	A high level of efficiency and a low level of inflammation.
System blend	Galvanic + High frequency	Combined action of electrolysis and thermolysis. A return electrode is required.	A high level of efficiency and a high level of inflammation.





The energy is channelled through a delicate needle, which is inserted through the follicle pore until it reaches the root where the electrical discharge is applied, causing irreversible damage to the hair germinal cells, thus eliminating the possibility of the growth of a new hair.


ELECTRICAL EPILATION À LA CARTE

Depilplus 13 MHz

Depilplus 13MHz is the electro-aesthetic equipment which provides the most effective electrical epilation techniques for the beauty professional. Prepares the skin and permanently removes unwanted hair in localised areas or hair resistant to photo-epilation.

Technical specifications

- | Depilplus 13MHz has pre-designed programmes which can also be personalised with different types of emission, intensities and time.
- | Emission of dry ozonized air in the needle holder with bactericidal and sedative effect.
- | Power level and timing regulator.
- | Pedal for controlling emission.



Technique Electrolysis

Thermolysis

Soft-Thermolysis

Current

Galvanic

Low power

Radio-frequency Thermal

TREATMENTS

Current

Chemical

 PERMANENT REMOVAL OF UNWANTED HAIR

Current

Very quick

Great effectiveness

Areas of medium sensitivity

	Soft-I nermolysis Flash	High power		Highly sensitive areas
	Blend	Galvanic+ Radio-frequency	Chemical + Thermal	Highly effective
	Thermal probe	Heat	Thermal	Product penetration
	Anaphoresis Prob	e Galvanic	Vasodilator effect	Facilitates channelling
	Cataphoresis post-depilación	Galvanic	Neutralising	Regulates skin pH Probe
	W			
The state of the s				
		NOSE		==
		EDI		EU
Contract of the Contract of th	III ENERGIA	0.0	II	





Range of emulsions and conductor gels enriched with natural active ingredients









enriched with draining and anti-cellulite active ingredients

Recommended for Kestos T2 and TEI System

Ref. 660070P · 1000ml



Thermoactive Emulsion. **Prevents cellulite**

Recommended for Rollaction

Ref. 660080P · 1000ml



Dry oil

Recommended for Rollaction

Ref. 660085P · 500ml

Our gels have passed all laboratory tests and guarantee care of the probes and energy applicators.

Our container, with dispenser included, is perfect for use in both facial and body treatments.

















Conductive gel for IPL and ultrasounds

Recommended for US1

Ref. 6600051 · 1000ml Ref. 6600050 · 5000ml



Conductive emulsion for RF

Recommended for Mesobiolift and My Sekret

Ref. 660095P · 1000ml



Conductive emulsion for endothermal wave. Firming active ingredients

Recommended for Mesobiolift RF

Ref. 660110P · 1000ml



Conductive emulsion for endothermal wave. Decongestant ingredients

Recommended for Mesobiolift RF

Ref. 660100P · 1000ml



Conductive emulsion for endothermal energy enriched with stem cells

Recommended for Mesobiolift RF

Ref. 660092P · 500ml



Conductive emulsion for capacitive and resistive treatment

Recommended for Binare Premium

Ref. 660200P · 500ml



Conductive gel for aesthetic electrostimulation

Recommended for Binare Care

Ref. 660210P · 1000ml



Rollaction System REF.0028460

DIMENSIONS	Width Height Depht	490 mm 950 mm 380 mm
WEIGHT BETA STAND REF. (VOLTAGE IR LIGHT MAX. SPEED	0077220 (Included)	19 kg 790x450 mm 115 W 850 nm 140 rpm

CODE	ACCESSORY Q	UANTIT
04RA011	ROLL Applicator 16mm	1
04RA012	ROLL Applicator 21mm	1
04RA013	ROLL Applicator 30mm	1
04RA014	ROLL Applicator 45mm	1
6600080	Thermoactive emulsion (1000m	nl) 1
6600085	Lipo Dren (500ml)	1

04EP013

Press G3		REF.0028500	-
DIMENSIONS		590 mm hout stand 180 mm 280 mm	166
WEIGHT BETA STAND POWER CONS	REF. 0077220 (Optiona	9 kg l) 790x450 mm 100 W	
CODE 0406700 0406710 0406711 0406720 0406721 0406730 0406731	ACCESSORY Abdomen Right leg Left leg Right foot Left foot Right arm Left arm Sleeve	QUANTITY 1 1 1 1 1 1 1 1	

Disposable trousers



REF.0080000



Kestos T2

DIMENSIONS

33TE114

REF.0028510

590 mm

	Height without star Depht	180 mm 280 mm
WEIGHT BETA STAND	REF. 0077220 (Optional) 79	10 kg 20x450 mm 300 W
CODE	ACCESSORY	QUANTITY
24TE235	Thermal bands 35W	2
24TE125	Thermal bands 25W	4
0400800	Elastic bands 100cm	4
0400900	Elastic bands 80cm	4
0401000	Elastic bands 60cm	4
0401100	Elastic bands 40cm	4
33TE114	Band cables	6
6600071	Conductive gel - Firming 1.00	10 ml 1
6600070	Conductive gel - Draining 1.00	10 ml 1
30T0110	Hanger bands	1
	COURTESY KIT	
24TE235	Thermal bands 35 W	1
24TE125	Thermal bands 25 W	1

Band cables

Width



US1 REF.0028490 DIMENSIONS Width 310 mm Height 350 mm Depht 180 mm WEIGHT 2,5 kg FREQUENCY 1MHz POWER 0,5 - 1,8W MODE Pulsed and continuous POWER CONSUMPTION 90W ACCESSORY QUANTITY CODE 0403992 1-C180 probe 0403998 2-C270 probe 6600051 Conductive gel 1.000ml

My Sekr	et	REF.0710000	1
DIMENSIONS	Width Height Depht	310 mm 350 mm 180 mm	
WEIGHT FREQUENCY POWER MODE POWER CONSUM	PTION	2,5 kg 3MHz 30W Continuous 90W	
40SK120 E 40SK110 F	ACCESSORY Body electrode Facial electrode RF Emulsion 1.000ml	QUANTITY 1 1	

TECHNICAL DATA

Binary Premium REF.0089020



TermoBeauty REF.0089040

Complet

REF.0028100

DIMENSIONS	Height without stand	490 mm 180 mm 280 mm	All
WEIGHT BETA STAND REF		5,75 kg 450 mm 69 mw	
	ACCESSORY HIGH FREQUENCY	QUANTITY	
3700700 M 3700500 F 3700600 F	Electrode grip Neck electrode Facial electrode Pencil electrode	1 1 1	
l	ndirect treatment electrode	1	E A De
4000700 N	Electrodes with grip Metal electrode Deep cleansing electrode	2 1 1	
0401320 E	Rollers Bags x 30 disposable electrodes Red – black cable	2 2 1	
	Multiple cable (5 red 5 black) Conductive gel (1000ml)	1 1	
9	SPRAY / SUCTION		
	Silicone tubes 5 x 8 x 1500 mm Round suction cap	2	
040270*	Spray bottles Flat suction cap	2	
E	BROSSAGE		
	Brush grip Nylon brush ø 20 mm	1 1	
	Nylon brush ø 40 mm Nylon brush ø 40 mm	1 1	

Depilplus 13 MHz REF.0028600

DIMENSIONS		490 mm It stand 180 mm 280 mm	Jwl
WEIGHT BETA STAND	REF. 0077220 (Optional)	5 kg 790x450 mm	SEPHEN
POWER CON FREQUENCY POWER	SUMPTION	450mA 50 Hz 105W	
CODE	ACCESSORY	QUANTITY	
0401700	DIN pedal	1	
040181*	Needle holde	1	
040182*	Thermal probe	1	
0407301	Straight needle 50 units	box. 1	
3306140	Black return cable	1	Ave.
4006909	Large stainless steel ret	urn electrode 1	LA Jan

Depilplus 3

REF.0028610

DIMENSIONS	Width Height Depht	310 mm 350 mm 180 mm
WEIGHT POWER CONSUMPTION FREQUENCY POWER		2,5 kg 450 mA 50 Hz 105 W

CODE	ACCESSORY	QUANTITY	
0401700	DIN pedal	1	
040181*	Needle holder	1	
0407301	Straight needle 50 units box	1	
3306140	Black return cable	1	
4006909	Large stainless steel return	electrod	1

OUR AESTHETIC TECHNIQUE'S GUIDE

RÖS'S provides health and beauty experts the most suitable technological tools for each treatment. Some of them have become, today, classics; others, are the last in technical innovations. Whatever they are the requirements you are looking for your beauty salon, we offer you a proper solution to achieve the best results.

Choose those with whom you find comfortable yourself!

BODY TECHNIQUES	ENERGY	APPLICATOR	MAIN BENEFITS	RÖS'S EQUIPMENT
THERMO-STIMULATION	Thermotherapy + electrostimulation	Thermo-stimulation strips	Reduction, sculpting and firming	TEI System, Kestos T2
ULTRASOUNDS	Ultrasounds	1MHz or 3MHz probe	Cellulite treatment	US1
SONOPHORESIS	Ultrasounds	1MHz or 3MHz probe	Reduction, sculpting and firming	US1
PHYSIO-ACTIVE MASSAGE	Compresion massage	Rollers	Reduction, sculpting and firming	Rollaction System
PRESSOMASSAGE	Pneumatic massage	Legs, abdomen and arms	Draining, activation of circulation	Press G3
ELECTROLYSIS / TERMOLYSIS	Electrical epilation	Needle	Hair removal	Depilplus 13MHz, Depilplus 3
RADIOFREQUENCY	Radiofrequency	Body Electrode	Firming	Mesobiolift RF, My Sekret
ENDOTHERMAL WAVE	RF + Direct current	Body Electrode	Firming and decongestant treatments	Mesobiolift RF
DIATHERMY	RF Capacitive and resistive	Body electrode	Cell regeneration	Binary Premium TermoBeauty
ELECTRO-STIMULATION	Currents	Stimulations strips	Firming and lipolysis	Binary Premium, TEI, Kestos T2

FACIAL TECHNIQUES	ENERGY	APPLICATOR	MAIN BENEFITS	RÖS'S EQUIPMENT
BRUSHES	Mechanic brushes	Brushes	Cleaning	Complet
HIGHT FREQUENCY	Alternating current	Electrode (direct/indirect)	Bactericide and antisept	Complet
DESCALING	Galvanic current	Rubber electrode	Deep cleansing	Complet
IONISATION	Galvanic current	Electrodes or rollers	Introduction of ionisable cosmetics	Complet, Mesobiolift RF
ELECTROLIFTING	Low frequency micro-currents	Electrodes	Muscle and tissues firming lifting effect	Complet, Binary Premium
ENDOTHERMAL WAVE	RF + Direct current	Facial electrode	Firming and restoring	Mesobiolift RF
DIATHERMY	RF Capacitive and resistive	Facial electrode	Cell regeneration	Binary Premium TermoBeauty



